



UNIPLAST

From Lab To Process Applications

Plastic Lab Ware Catalogue






From Lab To Process Applications



Celebrating our Past

Inspiring the Future



Welcome to the Latest Edition of UNIPLAST Plastic Lab ware catalog

Uniplast was established in 1969 in New Delhi, India, started our work to serve science and medical field. We began manufacturing and exporting science and medical products with high precision. At Uniplast, we have pioneered in General and special Medical / Scientific equipment, be it to Laboratories, Research Center, Medicals, Technical Institute, Universities, Schools, Colleges etc. In the year 1998, Uniplast had diversified into the manufacturing of PVC laminated sheets and PET manufacturing. This catalog only demonstrates a real depth in our product range related to Scientific product range and covers a large aspect of the end user's needs. We consistently keep evaluating and include new products into our Catalog, keeping in mind new areas of research and the need of all the end users.

Manufacturing capabilities: All our products are crafted using best grade raw materials and advanced technology in accordance with industry standards. At Uniplast, we are fully equipped with Knowledge base, latest machinery, clean room capability, infrastructure and tools to design and develop new custom products as per the need of our customers. This unit helps us in meeting the individual requirements of our clients in the best possible manner.

Mission: Our Mission has always been to offer high quality products to serve our wide range of scientific applications across sectors such as Health care, Education, Pharmaceuticals, Hospitals, Biotechnology, and Research & Academics.

We would like to thank our customers and channel partners for their continued support and hope that Uniplast products would continue to meet their expectations and be the first choice to fulfill their future needs.

Shiksha Tyagi

Director



Regd office: BP 174, Shalimar Bagh (West), Delhi 110088, INDIA

Works : D 221, Sector 3, DSIIIDC, Bawana Industrial Area, New Delhi -110039, INDIA

Tele: +91-9971700758, +91-9971700756, +91-9717867696, +91-11-27479180, +91-11-27762597

Email: uniplast@ymail.com, uniplast@uniplast.us

Website: <http://uniplast.us>

Index

Bottles

Aspirator Bottles	01
Dropper Bottles T.K, With Lid	01
Wash Bottles, Usual Type & With Dispensor	02
Self Venting Wash Bottles	02
Wash Bottles Wide Mouth	03
Narrow Mouth Bottles (PP, HDPE, AMBER, LDPE)	03
Wide Mouth Bottles (PP, HDPE, AMBER, LDPE)	04
Narrow Mouth Bottles Print Graduation	05
Wide Mouth Bottles Print Graduation	05
Wide Mouth Chemical Bottles	06
Narrow Mouth Chemical Bottles	06
Animal Feeder Water Bottles	06
Storage Bottles	06
Squeeze Bottle Dispensor	07
Multiuse Container	07

Volumetric's

Measuring Cylinder	08
Measuring Jug	08
Measuring Jug Printed	08
Measuring Beaker	09
Measuring Beaker Printed	09
Volumetric Flask	07
Conical Flask	07
Over Flow Vessel	19
Rain Gauze Cylinder With Collector	09

Micro Biology

Petri Dish	12
Coplin Jar	13
Slide Boxes	13
Specimen Jar	13
Desiccator	12
Micro Centrifuge Tube	14
Ria Vial	14
Cryo Vial	14
Centrifuge Tube Conical Bottom	14
Spatula	17
Stirring Rods, Inoculation Needle & loop	17
Cryo Box	16

Burette & Funnel

Burettes	10
Retort Stand	10
Burette Clamp	10
Funnel Holder	10
Separatory Funnel Holder	11
Separatory Funnel	11
Analytical Funnel	11

Stands & Racks

Conical Centrifuge Tubes Rack	15
Round Test Tube Stand	15
Rectangular Test Tube Stand	15
Micro Tubes Rack	15
Test Tube Stand 3 (Tier)	15

Containers

Pneumatic Trough	18
Simple Cell Pot	18
Leclanche Cell Pot	18
Laboratory Tray	18
Plantation Pots	19
Urine / Stool Container	19
Urine Pot (Male /Female)	19

School & Science Kits

Kipps Apparatus	11
Paschal Law Apparatus	20
Siphon Automatic Apparatus	20
Atomic Model Set	20
Crystal Model Set	20
Connectors	20

Pipette Products

Pipette Pump	16
Pasture Pipettes	16
Micro Pipettes	17
Pipette Stand Vertical	16
Pipette Stand Horizontal	16

PLASTIC LABWARE

1 ASPIRATOR BOTTLES

Material: Polypropylene (P.P.)
Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
101114	1.5 Ltr.	53 MM.	180 MM.	2 Pcs.	420.00
101115	05 Ltr.	53 MM.	357 MM.	1 Pcs.	600.00
101116	10 Ltr.	53 MM.	465 MM.	1 Pcs.	1070.00
101117	20 Ltr.	53 MM.	545 MM.	1 Pcs.	1520.00
101119	Spare stop cock	—	—	1 Pcs.	450.00



2 DROPPER BOTTLES (T. K. Type)

These Dropping Bottles are leak proof & Perfection to dispense uniform volume of liquid when passed lightly.

Material: Low Density Polyethylene (L.D.P.E.)
Color: Natural White, Translucent

Code	Capacity	I.D. Neck	Height	Packing	Price Per Pc. (Rs.)
101201	30 ML.	15.5 MM.	72 MM.	36 Pcs.	11.00
101202	60 ML.	16.0 MM.	91 MM.	36 Pcs.	13.00
101203	125 ML.	16.0 MM.	91 MM.	12 Pcs.	17.00



3 Dropping Bottles with Lid

Material: Low Density Polyethylene (L.D.P.E.)
Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
101305	30 ML.	22.5 MM.	75 MM.	36 Pcs.	15.00
101306	60 ML.	22.5 MM.	100 MM.	36 Pcs.	16.50
101307	125 ML.	22.5 MM.	112 MM.	12 Pcs.	19.50
101308	250 ML.	22.5 MM.	145 MM.	12 Pcs.	31.50



4 WASH BOTTLE

Material: Low Density Polyethylene (L.D.P.E.)
Color: Natural White, Translucent



Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
101401	125 ML.	24 MM.	100 MM.	12 Pcs.	18.00
101402	250 ML.	24 MM.	133 MM.	12 Pcs.	26.00
101403	500 ML.	28 MM.	170 MM.	12 Pcs.	35.00
101404	1000 ML.	38 MM.	201 MM.	06 Pcs.	60.00

5 WASH BOTTLE SIDE DISPENSOR

Material: Low Density Polyethylene (L.D.P.E.)
Color: Natural White, Translucent



Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
101501	125 ML.	24 MM.	120 MM.	12 Pcs.	26.50
101502	250 ML.	38 MM.	138 MM.	12 Pcs.	43.50
101503	500 ML.	38 MM.	170 MM.	12 Pcs.	51.00

6 SELF VENTING LABELED WASH BOTTLE

Material: Low Density Polyethylene (L.D.P.E.)
Color: Natural White, Translucent



Code	Capacity	Chemical	Screw O.D.	Packing	Price Per Pc. (Rs.)
101710	500 ML.	Acetone	53 MM.	06 Pcs.	175.00
101711	500 ML.	Distilled Water	53 MM.	06 Pcs.	175.00
101712	500 ML.	Isopropanol	53 MM.	06 Pcs.	175.00
101713	500 ML.	Ethanol	53 MM.	06 Pcs.	175.00
101714	500 ML.	Methanol	53 MM.	06 Pcs.	175.00
101715	500 ML.	Sodium Hypo Chlorite	53 MM.	06 Pcs.	175.00
101716	500 ML.	Assorted	53 MM.	06 Pcs.	175.00
101810	1000 ML.	Acetone	53 MM.	06 Pcs.	275.00
101811	1000 ML.	Distilled Water	53 MM.	06 Pcs.	275.00
101812	1000 ML.	Isopropanol	53 MM.	06 Pcs.	275.00
101813	1000 ML.	Ethanol	53 MM.	06 Pcs.	275.00
101814	1000 ML.	Methanol	53 MM.	06 Pcs.	275.00
101815	1000 ML.	Sodium Hypo Chlorite	53 MM.	06 Pcs.	275.00
101816	1000 ML.	Assorted	53 MM.	06 Pcs.	275.00

7 WASH BOTTLE WIDE MOUTH

Material: Low Density Polyethylene (L.D.P.E.)

Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
101601	500 ML.	53 MM.	170 MM.	06 Pcs.	88.00
101602	1000 ML.	53 MM.	200 MM.	06 Pcs.	150.00



8 NARROW MOUTH REAGENT BOTTLE - PP

Material: Polypropylene (P.P.), Autoclavable

Color: Natural, Transparent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102103	30 ML.	21 MM.	62 MM.	36 Pcs.	09.00
102104	60 ML.	21 MM.	85 MM.	24 Pcs.	11.00
102105	125 ML.	24 MM.	101 MM.	12 Pcs.	14.50
102106	250 ML.	24 MM.	133 MM.	12 Pcs.	24.50
102107	500 ML.	28 MM.	171 MM.	12 Pcs.	30.50
102108	1000 ML.	38 MM.	201 MM.	06 Pcs.	59.00



9 NARROW MOUTH REAGENT BOTTLE - HDPE

Material: High Density Polyethylene (H.D.P.E.)

Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102203	30 ML.	21 MM.	62 MM.	36 Pcs.	08.00
102204	60 ML.	21 MM.	85 MM.	24 Pcs.	09.50
102205	125 ML.	24 MM.	101 MM.	12 Pcs.	13.50
102206	250 ML.	24 MM.	133 MM.	12 Pcs.	22.00
102207	500 ML.	28 MM.	171 MM.	12 Pcs.	27.00
102208	1000 ML.	38 MM.	201 MM.	06 Pcs.	54.00



10 NARROW MOUTH REAGENT BOTTLE - AMBER

Material: High Density Polyethylene (H.D.P.E.)

Color: Amber, Opaque

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102303	30 ML.	21 MM.	62 MM.	36 Pcs.	09.50
102304	60 ML.	21 MM.	85 MM.	24 Pcs.	11.50
102305	125 ML.	24 MM.	101 MM.	12 Pcs.	16.00
102306	250 ML.	24 MM.	133 MM.	12 Pcs.	27.50
102307	500 ML.	28 MM.	171 MM.	12 Pcs.	33.50
102308	1000 ML.	38 MM.	201 MM.	06 Pcs.	54.00



11 NARROW MOUTH REAGENT BOTTLES - LDPE

Material : Low Density Polyethylene (L.D.P.E.)
Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102403	30 ML.	21 MM.	62 MM.	36 Pcs.	08.50
102404	60 ML.	21 MM.	85 MM.	24 Pcs.	10.50
102405	125 ML.	24 MM.	101 MM.	12 Pcs.	13.50
102406	250 ML.	24 MM.	133 MM.	12 Pcs.	23.50
102407	500 ML.	28 MM.	171 MM.	12 Pcs.	29.00
102408	1000 ML.	38 MM.	201 MM.	06 Pcs.	56.00



12 WIDE MOUTH REAGENT BOTTLE - PP

Material : Polypropylene, Autoclavable
Color: Natural, Transparent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102603	30 ML.	28 MM.	62 MM.	36 Pcs.	10.00
102604	60 ML.	28 MM.	85 MM.	24 Pcs.	11.50
102605	125 ML.	38 MM.	100 MM.	12 Pcs.	16.50
102606	250 ML.	43 MM.	130 MM.	12 Pcs.	28.50
102607	500 ML.	53 MM.	170 MM.	12 Pcs.	36.50
102608	1000 ML.	63 MM.	199 MM.	06 Pcs.	65.00



13 WIDE MOUTH REAGENT BOTTLES - HDPE

Material : High Density Polyethylene (H.D.P.E.)
Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102703	30 ML.	28 MM.	62 MM.	36 Pcs.	09.50
102704	60 ML.	28 MM.	85 MM.	24 Pcs.	11.50
102705	125 ML.	38 MM.	100 MM.	12 Pcs.	14.50
102706	250 ML.	43 MM.	130 MM.	12 Pcs.	25.50
102707	500 ML.	53 MM.	170 MM.	12 Pcs.	33.50
102708	1000 ML.	63 MM.	199 MM.	06 Pcs.	58.00



14 WIDE MOUTH REAGENT BOTTLES - AMBER

Material : High Density Polyethylene (H.D.P.E.)
Color: Amber, Opaque

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102803	30 ML.	28 MM.	62 MM.	36 Pcs.	11.50
102804	60 ML.	28 MM.	85 MM.	24 Pcs.	13.50
102805	125 ML.	38 MM.	100 MM.	12 Pcs.	18.00
102806	250 ML.	43 MM.	130 MM.	12 Pcs.	32.00
102807	500 ML.	53 MM.	170 MM.	12 Pcs.	40.00
102808	1000 ML.	63 MM.	199 MM.	06 Pcs.	62.00

15 WIDE MOUTH REAGENT BOTTLES - LDPE

Material : Low Density Polyethylene (L.D.P.E.)
Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
102903	30 ML.	28 MM.	62 MM.	36 Pcs.	09.50
102904	60 ML.	28 MM.	85 MM.	24 Pcs.	11.50
102905	125 ML.	38 MM.	100 MM.	12 Pcs.	15.50
102906	250 ML.	43 MM.	130 MM.	12 Pcs.	27.50
102907	500 ML.	53 MM.	170 MM.	12 Pcs.	34.50
102908	1000 ML.	63 MM.	199 MM.	06 Pcs.	62.00

16 NARROW MOUTH BOTTLE PRINT GRADUATION

Material : Polypropylene (PP), Autoclavable
Color: Natural, Transparent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
103501	125 ML.	24 MM.	101 MM.	12 Pcs.	27.00
103502	250 ML.	24 MM.	133 MM.	12 Pcs.	36.00
103503	500 ML.	28 MM.	171 MM.	12 Pcs.	45.00
103504	1000 ML.	38 MM.	201 MM.	06 Pcs.	76.00

17 WIDE MOUTH BOTTLE PRINT GRADUATION

Material : Polypropylene (PP), Autoclavable
Color: Natural, Transparent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
103801	125 ML.	38 MM.	100 MM.	12 Pcs.	28.00
103802	250 ML.	43 MM.	130 MM.	12 Pcs.	37.00
103803	500 ML.	53 MM.	170 MM.	12 Pcs.	46.00
103804	1000 ML.	63 MM.	199 MM.	06 Pcs.	78.00



18 WIDE MOUTH CHEMICAL BOTTLE - AMIDA

These Bottles are suitable to store Insecticides & Pesticides Chemicals

Material : High Density Polyethylene (H.D.P.E.)

Color: White, Translucent



Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
103601	50 ML.	28 MM.	85 MM.	36 Pcs.	11 .00
103602	100 ML.	38 MM.	100 MM.	12 Pcs.	15 .00
103603	250 ML.	43 MM.	130 MM.	12 Pcs.	25 .00
103604	500 ML.	53 MM.	170 MM.	12 Pcs.	33 .00
103605	1000 ML.	63 MM.	220 MM.	06 Pcs.	60 .00

19 NARROW MOUTH CHEMICAL BOTTLE - AMIDA

These Bottles are suitable to store Insecticides & Pesticides Chemicals

Material : High Density Polyethylene (H.D.P.E.)

Color: White, Translucent



Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
103701	50 ML.	21 MM.	85 MM.	36 Pcs.	11 .00
103702	100 ML.	24 MM.	101 MM.	12 Pcs.	15 .00
103703	250 ML.	24 MM.	133 MM.	12 Pcs.	25 .00
103704	500 ML.	28 MM.	171 MM.	12 Pcs.	33 .00
103705	1000 ML.	38 MM.	220 MM.	06 Pcs.	60 .00

20 WATER BOTTLE ANIMAL FEEDER

Material : Polypropylene (PP), Autoclavable

Color: Natural, Transparent



Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
103901	125 ML.	24 MM.	170 MM.	12 Pcs.	64 .00
103902	250 ML.	24 MM.	210 MM.	12 Pcs.	94 .00

21 STORAGE BOTTLES

Material : Polypropylene (PP), Autoclavable

Color: Natural, Transparent



Code	Capacity	Shape	Screw O.D.	Packing	Price Per Pc. (Rs.)
103301	250 ML.	Rectangular	24 MM.	12 Pcs.	10 .00
103302	1000 ML.	Rectangular	38 MM.	06 Pcs.	56 .00
103303	200 ML.	Round	65 MM.	12 Pcs.	12 .00

22 CONICAL FLASK

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
103001	125 ML.	26 MM.	128 MM.	12 Pcs.	34.00
103002	250 ML.	38 MM.	150 MM.	12 Pcs.	50.00
103003	500 ML.	38 MM.	173 MM.	12 Pcs.	80.00



23 VOLUMETRIC FLASK

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
104101	25 ML.	15.5 MM.	136 MM.	12 Pcs.	44.00
104102	50 ML.	16.5 MM.	144 MM.	06 Pcs.	61.00
104103	100 ML.	18.5 MM.	160 MM.	06 Pcs.	74.00
104104	250 ML.	19.0 MM.	225 MM.	06 Pcs.	92.00
104105	500 ML.	20.0 MM.	270 MM.	03 Pcs.	131.00
104106	1000 ML.	21.0 MM.	310 MM.	03 Pcs.	196.00



24 SQUEEZE BOTTLE DISPENSOR

Material : Low Density Polyethylene (L.D.P.E)
Color: Natural White, Translucent

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
104001	125 ML.	35 MM.	170 MM.	12 Pcs.	22.00
104002	250 ML.	35 MM.	184 MM.	12 Pcs.	26.00
104003	500 ML.	37 MM.	188 MM.	12 Pcs.	34.00



25 MULTI USE CONTAINER

Material : Poly Vinyl Chloride (P.V.C.)
Color: Natural Clear

Code	Capacity	Screw O.D.	Height	Packing	Price Per Pc. (Rs.)
104008	1500 ML.	50 MM.	190 MM.	06 Pcs.	50.00



26 MEASURING CYLINDER

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent



Code	Capacity	Least Count	Packing	Price Per Pc. (Rs.)
113601	10 ML.	01	24 Pcs.	11.50
113602	25 ML.	01	12 Pcs.	18.00
113603	50 ML.	01	12 Pcs.	22.00
113604	100 ML.	01	12 Pcs.	26.50
113605	250 ML.	05	12 Pcs.	55.00
113606	500 ML.	05	12 Pcs.	84.00
113607	1000 ML.	05	06 Pcs.	170.00
113608	2000 ML.	50	01 Pcs.	440.00

27 MEASURING JUG

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent



Code	Capacity	Least Count	Packing	Price Per Pc. (Rs.)
113902	250 ML.	05	12 Pcs.	22.00
113903	500 ML.	20	12 Pcs.	27.00
113904	1000 ML.	100	06 Pcs.	48.00
113905	2000 ML.	100	06 Pcs.	85.00

28 MEASURING JUG PRINTED

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent



Code	Capacity	Least Count	Packing	Price Per Pc. (Rs.)
114002	250 ML.	50	12 Pcs.	34.50
114003	600 ML.	100	12 Pcs.	54.00
114004	1000 ML.	100	06 Pcs.	84.00

29 MEASURING BEAKER

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent

Code	Capacity	Least Count	Packing	Price Per Pc. (Rs.)
114201	50 ML.	05	24 Pcs.	07.50
114202	100 ML.	05	12 Pcs.	08.50
114203	250 ML.	05	12 Pcs.	17.00
114204	600 ML.	20	12 Pcs.	26.00
114205	1000 ML.	100	06 Pcs.	44.00
114206	2000 ML.	100	06 Pcs.	87.00
114207	5000 ML.	1000	03 Pcs.	250.00



30 MEASURING BEAKER PRINTED

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent

Code	Capacity	Least Count	Packing	Price Per Pc. (Rs.)
114301	50 ML.	10	24 Pcs.	14.00
114302	100 ML.	20	12 Pcs.	15.50
114303	250 ML.	50	12 Pcs.	27.00
114304	600 ML.	100	12 Pcs.	41.50
114305	1000 ML.	100	06 Pcs.	80.00



31 RAIN GAUZE CYLINDER

These Cylinder are available in Transparent and Translucent with size 200 cm square useful in Gauzing Rain water.

Material Option : Polypropylene (Translucent) or Polystyrene (Transparent)

Code	Material	Capacity	Packing	Price Per Pc. (Rs.)
114501	P.P.	200 ml.	12 Pcs.	38.00
114502	P.S.	200 ml.	12 Pcs.	42.00

32 RAIN GAUZE CYLINDER WITH COLLECTOR

Complete unit is Transparent , Cylinder is made of Poly Styrene and Collector is made of clear PVC.

Code	Material	Capacity	Packing	Price Per Pc. (Rs.)
114503	P.S. / P.V.C.	200 ml.	01 Pc.	200.00






33 RETORT STAND

Material : A.B.S. (Acrylonitrile Butadiene Styrene)
Color: White


Code	Description	Packing	Price Per Pc. (Rs.)
125401	225 x 150 mm, Side Hole	01 Pc.	335.00
125402	300 x 225 mm, Center Hole	01 Pc.	560.00



34 BURETTE

Material : Polymethylpentene TPX

Code	Capacity	Packing	Price Per Pc. (Rs.)
125801	25 ML.	10 Pcs.	230.00
125802	50 ML.	10 Pcs.	250.00



35 BURETTE CLAMP

Material : Polypropylene (P.P.)
Color: White

Code	Description	Packing	Price Per Pc. (Rs.)
124501	Burette Clamp Single	12 Pcs.	50.00
124502	Burette Clamp Double	12 Pcs.	88.00



36 FUNNEL HOLDER

Material: Polypropylene (P.P.)

Code	Description	Packing	Price Per Pc. (Rs.)
124504	Funnel Holder Single	12 Pcs.	48.00

37 SEPARATORY FUNNEL HOLDER

Material : Polypropylene (P.P.)
Color: White

Code	Description	Packing	Price Per Pc. (Rs.)
124507	Single	12 Pcs.	54.00



38 SEPARATORY FUNNEL

Material : Polypropylene (P.P.), Autoclavable

Code	Capacity	Packing	Price Per Pc. (Rs.)
145101	100 ML.	2 Pcs.	164.00
145102	250 ML.	2 Pcs.	194.00
145103	500 ML.	2 Pcs.	240.00



39 ANALYTICAL FUNNELS

Material : Polypropylene (P.P.), Autoclavable

Code	Particulars	Inside	Stem Length	Packing	Price Per Pc. (Rs.)
145302	50 MM.	49MM.	39	36 Pcs.	06.00
145303	62 MM.	61 MM.	60	36 Pcs.	08.50
145304	75 MM.	74 MM.	70	36 Pcs.	09.50
145305	100 MM.	93 MM.	84	24 Pcs.	16.00



40 KIPP'S APPARATUS

Material: High Density Polyethylene (H.D.P.E.)

Code	Capacity	Packing	Price Per Pc. (Rs.)
156101	250 ML.	01 Pc.	575.00
156102	500 ML.	01 Pc.	660.00
156103	1000 ML.	01 Pc.	820.00
156104	Medium	01 Pc.	400.00
156105	Large	01 Pc.	500.00
Spare	Thistle Funnel Set For Kipp's	01 Pc.	50.00



41 DESICCATOR (PLAIN)

Material: Polypropylene / Polycarbonate



Code	Capacity	Packing	Price Per Pc. (Rs.)
205601	150 MM.	01 Pc.	900.00
205602	200 MM.	01 Pc.	1500.00
205603	250 MM.	01 Pc.	2100.00
205604	300 MM.	01 Pc.	3400.00
Spare	Silicon Ring	01 Pc.	450.00

42 DESICCATOR (VACUUM)

Material: Polypropylene / Polycarbonate



Code	Capacity	Packing	Price Per Pc. (Rs.)
205801	150 MM.	01 Pc.	1400.00
205802	200 MM.	01 Pc.	2100.00
205803	250 MM.	01 Pc.	2600.00
205804	300 MM.	01 Pc.	4000.00
Spare	Stop Cork	01 Pc.	450.00

43 PETRI DISH REUSABLE

Material : Polypropylene (P.P.), Autoclavable
Color: Natural, Transparent



Code	Description	Packing	Price Per Pc. (Rs.)
167101	50 MM.	36 Pcs.	07.50
167102	75 MM.	36 Pcs.	10.00
167103	100 MM.	36 Pcs.	12.00

44 PETRI DISH DISPOSABLE

Material : Polystyrene (P.S.)
Color: Natural, Clear



Code	Description	Sterile	Packing	Price Per Pc. (Rs.)
167106	90 MM.	No	36 Pcs.	08.00
167107	90 MM.	Yes	36 Pcs.	09.50

Individually Wrapped

45 COPLIN JAR


Material: Polypropylene (P.P.)

Code	Capacity	Packing	Price Per Pc. (Rs.)
167201	10 Slides	12 Pcs.	24.00


46 SLIDE BOX (DETACHABLE)

Material: Polystyrene (P.S.)

Code	Capacity	Packing	Price Per Pc. (Rs.)
167205	5 Slides	24 Pcs.	35.00
167206	25 Slides	24 Pcs.	51.00
167207	50 Slides	12 Pcs.	97.00
167208	100 Slides	06 Pcs.	150.00


47 SLIDE BOX (HINGED)

Material: Polypropylene (P.P.)

Code	Capacity	Packing	Price Per Pc. (Rs.)
167211	25 Slides	24 Pcs.	20.00
167212	50 Slides	12 Pcs.	42.00
167213	100 Slides	12 Pcs.	100.00


48 SPECIMEN JAR

Material: Poly Styrene (P.S.)

Code	Description	Packing	Price Per Pc. (Rs.)
167301	Cylindrical Shape	12 Pcs.	38.00





49 MICRO CENTRIFUGE TUBE

Material: Polypropylene (P.P)

Code	Capacity	Packing	Price Per Pc. (Rs.)
167420	0.5 ML.	1000 Pcs.	00.75
167421	1.5 ML.	500 Pcs.	00.90



50 RIA VIAL

Material : Polystyrene (P.S.)


Code	Capacity	Packing	Price Per Pc. (Rs.)
167408	12 x 75 MM.	100 Pcs.	01.15



51 CRYO VIAL

Material: Polypropylene (P.P)

Code	Capacity	Packing	Price Per Pc. (Rs.)
167411	1.8 ML.	100 Pcs.	04.50
167412	4.5 ML.	100 Pcs.	08.00



52 CENTRIFUGE TUBES CONICAL BOTTOM


Material: Polypropylene (P.P)

Code	Capacity	Sterile	Packing	Price Per Pc. (Rs.)
167401	15 ML.	No	100 Pcs.	08.50
167402	15 ML.	Yes	100 Pcs.	10.00
167403	50 ML.	No	50 Pcs.	10.50
167404	50 ML.	Yes	50 Pcs.	12.00

53 CONICAL CENTRIFUGE TUBE RACK

Material : Polystyrene (P.S.), Autoclavable


Code	Description	Packing	Price Per Pc. (Rs.)
124008	Wells for 15 ML. Tubes measures 17 MM Dia Wells for 50 ML. Tubes measures 28 MM Dia	02 Pcs.	575.00



54 TEST TUBE STAND (ROUND / RECTANGULAR)

Material : Polypropylene (P.P.), Autoclavable

Code	Description	Packing	Price Per Pc. (Rs.)
124001	Round St.- 19 MM x 25 MM Tubes	12 Pcs.	107.00
124002	Rectangular St.- 16 MM x 25 MM Tubes	12 Pcs.	24.00



55 MICRO TUBE RACK

Material : Polypropylene (P.P.), Autoclavable


Code	Description	Packing	Price Per Pc. (Rs.)
124003	12 MM x 20 MM Tubes	12 Pcs.	40.00



56 TEST TUBE STAND 3 (TIER)

Material : Polypropylene (P.P.), Autoclavable

Code	Description	Packing	Price Per Pc. (Rs.)
124101	13 MM x 31 Tubes	02 Pcs.	94.00
124103	16 MM x 31 Tubes	02 Pcs.	94.00
124105	18 MM x 12 Tubes	04 Pcs.	60.00
124106	20 MM x 20 Tubes	02 Pcs.	134.00
124107	25 MM x 18 Tubes	04 Pcs.	134.00






57 PIPETTE STAND (VERTICAL)

Material : Polypropylene (P.P.), Autoclavable

Code	Description	Packing	Price Per Pc. (Rs.)
124201	For 28 Pipettes, Dia 14 MM, Heavy Base	01 Pcs.	335.00
124202	For 28 Pipettes, Dia 14 MM, Light Base	06 Pcs.	240.00



58 PIPETTE STAND (HORIZONTAL)

Material : Polypropylene (P.P.), Autoclavable

Code	Description	Packing	Price Per Pc. (Rs.)
124203	For 12 Pipettes	06 Pcs.	95.00



59 CRYO BOX

These Boxes are suitable for 1.5 and 2.0 ml. Micro centrifuge tubes, Cryo Vials and tubes.

Material : Polypropylene (P.P.), Autoclavable

Code	Capacity	Packing	Price Per Pc. (Rs.)
124210	100 Wells	01 Pcs.	600.00



60 PIPETTE PUMP

Material : Polypropylene (P.P.) / Silicon, Autoclavable

Code	Capacity	Packing	Price Per Pc. (Rs.)
205201	10 ML.	01 Pc.	230.00
205202	25 ML.	01 Pc.	260.00



61 PASTURE PIPETTE

Material : Low Density Polyethylene (L.D.P.E.)

Code	Capacity	Packing	Price Per Pc. (Rs.)
205205	03 ML.	500 Pcs.	01.40

62 MICROPIPETTE-SINGLE CHANNEL-FIXED VOLUME

Material : PVDF / ABS, Autoclavable

Code	Description	Volume Range (µL)	Price Per Pc. (Rs.)
205210		05 µL	1400.00
205211	Fix Volume Digital	10 µL	1400.00
205212	Micro-pipette	100 µL	1400.00
205213		250 µL	1400.00
205214		400 µL	1400.00
205215		500 µL	1400.00
205216		1000 µL	1400.00

63 MICROPIPETTE-SINGLE CHANNEL-VARIABLE VOLUME

Material : PVDF / ABS, Autoclavable

Code	Description	Volume Range (µL)	Increment	Price Per Pc. (Rs.)
205221		0.5 - 10 µL	0.02	4000.00
205222		05 - 50 µL	0.1	4000.00
205223	Variable Volume Digital	10 - 100 µL	0.1	4000.00
205224	Micro-pipette with Tip	20 - 200 µL	0.2	4000.00
205225	Injector.	100 - 1000 µL	1.0	4000.00
205226		1000 - 5000 µL	10.0	4000.00
205227		1000 - 10000 µL	20.0	4000.00

64 SPATULA


Material: Polypropylene (P.P.)


Code	Description	Packing	Price Per Pc. (Rs.)
227802	150 MM.	12 Pcs.	07.00
227803	200 MM.	12 Pcs.	09.00

65 STIRRING RODS, INOCULATION-NEEDLE & LOOP

Material : Polypropylene (P.P.), Polystyrene - Autoclavable

Code	Description	Packing	Price Per Pc. (Rs.)
227806	0-10 MM x H-250 MM	12 Pcs.	14.00
227807	0-06 MM x H-245 MM	12 Pcs.	14.00
227808	Inoculation Needle	100 Pcs.	04.50
227809	Soft Loop - 10µl	100 Pcs.	04.50
227810	L-Spreader	50 Pcs.	13.50





66 PNEUMATIC TROUGH

Material : Polypropylene (P.P.), Autoclavable

Code	Description	Packing	Price Per Pc. (Rs.)
228001	0-200 MM X H-100 MM	06 Pcs.	87.00



67 SIMPLE CELL POT

Material: Polypropylene (P.P.)


Code	Description	Packing	Price Per Pc. (Rs.)
228008	Simple Cell Pot	12 Pcs.	44.00



68 LECLANCHE CELL POT

Material: Polypropylene (P.P.)

Code	Description	Packing	Price Per Pc. (Rs.)
228009	Leclanche Cell Pot	12 Pcs.	54.00



69 LABORATORY TRAY

Material : Polypropylene (P.P.), Autoclavable

Code	Description	Packing	Price Per Pc. (Rs.)
228003	250 MM X 300 MM X 50 MM	05 Pcs.	190.00
228004	415 MM X 300 MM X 75 MM	05 Pcs.	240.00
228005	430 MM X 350 MM X 75 MM	05 Pcs.	320.00

70 PLANTATION POTS

Material : Polypropylene (P.P.)

Code	Description	Packing	Price Per Pc. (Rs.)
227901	100 MM.	20 Pcs.	24.00
227902	150 MM.	20 Pcs.	42.00
227903	200 MM.	20 Pcs.	75.00



71 URINE / STOOL CONTAINER

Material: Polypropylene (P.P.)

Code	Description	Packing	Price Per Pc. (Rs.)
238201	Urine 30 ML.	100 Pcs.	05.00
238202	Urine 60 ML.	100 Pcs.	06.00
238203	Urine 100 ML.	100 Pcs.	09.00
238204	Stool 60 ML.	100 Pcs.	13.00
238205	Stool 30 ML.	100 Pcs.	12.00



72 URINE POT

Material: Polypropylene (P.P.)

Code	Description	Packing	Price Per Pc. (Rs.)
103401	Male Shape	50 Pcs.	30.00
103402	Female Shape	50 Pcs.	30.00




73 OVER FLOW VESSEL

Cylindrical Shaped, with an Over flow Spout used to Measure Liquids.

Material: Polypropylene (P.P.)

Code	Description	Packing	Price Per Pc. (Rs.)
113201	Cylindrical Shape	12 Pcs.	52.00






74 ATOMIC MODEL SET

Material: Polypropylene (P.P.)


Code	Description	Packing	Price Per Pc. (Rs.)
158701	60 Balls Junior Set	01 Pcs.	380.00
158702	120 Balls Medium Set	01 Pcs.	600.00
158703	360 Balls Senior Set	01 Pcs.	1400.00



75 CRYSTAL MODEL SET

Material: Polyethylene (P.E.)

Code	Description	Packing	Price Per Pc. (Rs.)
158804	Diamond	01 Pcs.	200.00
158805	Sodium Chloride	01 Pcs.	200.00
158806	Graphite	01 Pcs.	200.00



76 CONNECTORS

Material: Polypropylene (P.P.)

Code	Description	Packing	Price Per Pc. (Rs.)
227101	06 MM. T-Connector	36 Pcs.	06.00
227102	06 MM. Y-Connector	36 Pcs.	06.00
227103	06 MM. L-Connector	36 Pcs.	06.00



77 PASCHAL LAW APPARATUS

Material : Polypropylene (P.P.) & Polyethylene (P.E.)

Code	Description	Packing	Price Per Pc. (Rs.)
158907	Piston Cylinder Fitted With Ethylene Bulb	12 Pcs.	200.00



78 SIPHON AUTOMATIC

This Siphon is Leak Proof available stop cock used to transfer liquids. Acid and Alkaline Proof

Material : Polypropylene (P.P.) & Polyethylene (P.E.)

Code	Description	Packing	Price Per Pc. (Rs.)
158908	Acid & Alkaline Proof	12 Pcs.	150.00

CHEMICAL RESISTANCE FOR UNIPLAST PRODUCTS

CLASSES OF SUBSTANCES AT ROOM TEMPERATURE	HDPE	LDPE	PC	PP	PS
ACIDS, DILUTE OR WEAK	Green	Green	Green	Green	Green
ACIDS, STRONG AND CONCENTRATED	Green	Green	Orange	Green	Yellow
ALCOHOLS, ALIPHATIC	Green	Green	Blue	Green	Green
ALDEHYDES	Blue	Blue	Yellow	Blue	Orange
BASES	Green	Green	Orange	Green	Green
ESTERS	Blue	Blue	Orange	Blue	Orange
HYDROCARBONS, ALIPHATIC	Blue	Yellow	Yellow	Blue	Orange
HYDROCARBONS, AROMATIC	Blue	Yellow	Orange	Yellow	Orange
HYDROCARBOS, HALOGENATED	Yellow	Orange	Orange	Yellow	Orange
KETONES	Blue	Blue	Orange	Blue	Orange
OXIDIZING AGENTS, STRONG	Yellow	Yellow	Orange	Yellow	Orange

- 30 Days of constant exposure cause no damage. Plastics may even tolerate for years.
- Little or no damages after 30 days of constant exposure to the reagents.
- Some effect after after 7 days of constant exposure to the reagent like crazing, cracking, loss of strength or discoloration.
- Not recommended for continuous use. Immediate damage may occur.

PHYSICAL PROPERTIES OF RAW MATERIAL USED FOR UNIPLAST PRODUCTS

Resin	Max use Temp°C	Brittle-ness Temp°C	Transparency							Flexi-bility	Water Absorp-tion%
				Autocla-ving	Gas	Dry Heat	Radia-tion	Disin-fectant	Specific Gravity		
HDPE	120	-100	Translucent	No	Yes	No	Yes	Yes	0.95	Rigid	20.01
LDPE	80	-100	Translucent	No	Yes	No	Yes	Yes	0.92	Excel	<0.01
PC	135	-135	Clear	Yes+	Yes	No	No	Yes	1.20	Rigid	0.35
PP	135	0	Translucent	Yes	Yes	No	No	Yes	0.90	Rigid	<0.02
PS	90	20	Clear	No	Yes	No	Yes	Some	1.05	Rigid	0.05

Sterilization :

- Autoclaving (1210C, 15 psig for 20 minutes) Clean and rinse item with distilled water before autoclaving. Certain chemicals which have no appreciable effect on resins at room temperature may cause deterioration at autoclaving temperature.
 - Always completely disengage threads before autoclaving.
 - Gas - Ethylene oxide formaldehyde.
 - Dry heat -1600C for 120 minutes.
 - Disinfectant - Benzalkonium chloride, formalin, ethanol etc.
 - Radiation - gamma irradiated at 2.5 Mard with unstabilized plastic.
- Sterilizing reduces mechanical strength, Do not use PC vessels for vaccum application if they have been autoclaved

TERMS & CONDITIONS OF SALE

- 1 This Catalog Includes all “Uniplast” Brand Products.
- 2 GST Will be charged extra as applicable.
3. The Price for items in bulk requirement will be Quoted on request.
4. All Prices are FOR NEW Delhi.
5. Uniplast Responsibility ceases once the cargo is delivered to the carrier / transporter
6. Goods will be dispatched in standard packing or multiple thereof.
7. Capacity measures are for laboratory use only and not for commercial practice.
8. Uniplast reserves the right to change / modify any product without prior notice.
9. All disputes shall be settled with in the jurisdiction of New Delhi.
- 10.This Catalog is effective 1st April, 2018, & cancels all over previous Catalogs.



UNIPLAST

From Lab To Process Applications

Regd office: BP 174, Shalimar Bagh (West), Delhi 110088, INDIA

Manufacturing : D 221, Sector 3, DSIIIDC, Bawana Industrial Area, New Delhi -110039, INDIA

Tele: +91-9971700758, +91-9971700756, +91-9717867696, +91-11-27479180

Email : uniplast@ymail.com

Website: <http://uniplast.us>